


Tucker Sling™ Positioning Aid

Support for babies who need elevation.

- One-piece support system for babies who need to be in elevated positions.
- May be used for babies with GER/acid reflux or respiratory problems.
- Attaches to the mattress or wedge like a fitted sheet.
- Baby is held firmly in place with adjustable, easy-to-release straps.
- Keeps babies both upright and secure and prevents them from sliding to the foot of an elevated mattress.
- Item descriptions indicate size of mattress and weight of baby.
- Optional foam wedge can be used under a crib mattress to assist in elevating the baby.
- Single-patient use, disposable.
- Item descriptions indicate size of mattress and weight of baby.

The Tucker Sling is the exact same product I introduced at Children's Medical Ventures. SandBox Medical became the distributor in January 2016 to ensure these products remain available to customers. Tucker Sling is manufactured by the inventor, Terry Jarret.

Read McCarty, Co-Founder, SandBox Medical
 Founder, Past CEO Children's Medical Ventures



| Item # | Compares to Philips Item# | Description | Qty. Per Box |
|------------|---------------------------|------------------------------------|--------------|
| TS16030112 | 91022-A1 | Incubator; Infant <2 lbs. Flannel | 12 |
| TS16030212 | 91022-A2 | Incubator; Infant 2-5 lbs. Flannel | 12 |
| TS16030312 | 91022-B1 | Bassinet; Infant 3-7 lbs. Flannel | 12 |
| TS16030412 | 91022-A3 | Crib; Infant 2-5 lbs. Flannel | 12 |
| TS16030512 | 91022-B2 | Crib; Infant 3-7 lbs. Flannel | 12 |
| TS16030612 | 91022-C1 | Crib; Infant 6-16 lbs. Flannel | 12 |
| TS16030712 | 91022-C2 | Crib; Infant 6-16 lbs. Quilted | 12 |
| TS16030812 | 91022-D | Crib; Infant 16-24 lbs. Quilted | 12 |
| TS16030912 | 91022-E | Crib; Infant 24-30 lbs. Quilted | 12 |
| TS1603101 | 91022-W | Foam wedge 12" x 26" x 24" | 1 |

Specifications

Sling: 100 % cotton flannel.

Wash in warm or cold water. Tumble dry, low heat.

To prevent lint from collecting in Velcro, fasten strips together before washing and drying.

Foam Wedge: Soy-based. Cover is nylon institutional mattress ticking; flame retardant, bacteria and stain resistant, anti-static, fluid proof.

Manufactured in the USA.

Philips is a registered trademark of Koninklijke Philips N.V.

5-1603-001B