


Tucker Sling™ Positioning Aid

Support for babies who need elevation.

- One-piece support system for babies who need to be in elevated positions.
- May be used for babies with GER/acid reflux or respiratory problems.
- Attaches to the mattress or wedge like a fitted sheet.
- Baby is held firmly in place with adjustable, easy-to-release straps.
- Keeps babies both upright and secure and prevents them from sliding to the foot of an elevated mattress.
- Item descriptions indicate size of mattress and weight of baby.
- Optional foam wedge can be used under a crib mattress to assist in elevating the baby.
- Single-patient use, disposable.
- Item descriptions indicate size of mattress and weight of baby.

The NOVAPLUS Tucker Sling is the exact same product I introduced at Children's Medical Ventures. SandBox Medical became the distributor in January 2016 to ensure these products remain available to Vizient members. Tucker Sling is manufactured by the inventor, Terry Jarret.

Read McCarty, Co-Founder, SandBox Medical
 Founder, Past CEO Children's Medical Ventures



Item #	Compares to Philips Item#	Description	Qty. Per Box
VPTS16030112	91022-A1	Incubator; Infant <2 lbs. Flannel	12
VPTS16030212	91022-A2	Incubator; Infant 2-5 lbs. Flannel	12
VPTS16030312	91022-B1	Bassinet; Infant 3-7 lbs. Flannel	12
VPTS16030412	91022-A3	Crib; Infant 2-5 lbs. Flannel	12
VPTS16030512	91022-B2	Crib; Infant 3-7 lbs. Flannel	12
VPTS16030612	91022-C1	Crib; Infant 6-16 lbs. Flannel	12
VPTS16030712	91022-C2	Crib; Infant 6-16 lbs. Quilted	12
VPTS16030812	91022-D	Crib; Infant 16-24 lbs. Quilted	12
VPTS16030912	91022-E	Crib; Infant 24-30 lbs. Quilted	12
VPTS1603101	91022-W	Foam wedge 12" x 26" x 24"	1

Specifications

Sling: 100 % cotton flannel.

Wash in warm or cold water. Tumble dry, low heat.

To prevent lint from collecting in Velcro, fasten strips together before washing and drying.

Foam Wedge: Soy-based. Cover is nylon institutional mattress ticking; flame retardant, bacteria and stain resistant, anti-static, fluid proof.

Manufactured in the USA.

Philips is a registered trademark of Koninklijke Philips N.V.

NP5-1603-001b

vizient™
 Awarded Supplier

2014 Pediatric Award Winner

NOVAPLUS®

NOVAPLUS is a registered trademark of Vizient Inc.

Contract number MS3961