

## 7 x 7 Inch Star-Shaped Pacifier Holder with Removable Premie or Newborn Size JollyPop® Pacifier

### PRODUCT SPECIFICATIONS:

Phthalate-free, latex-free, BPA-free and DEHP-free. Manufactured in the USA.

**CuddleStar:** Polyester non-woven outer and inner layer. Velcro® closure for pacifier. New, polypropylene beads contained in inner layer. Meets ASTM F963-11 safety code requirements and is baby safe.

**JollyPop Pacifier:** 100% medical grade silicone. Complies with USA Standard 16CFR 1511.

### INSTRUCTIONS FOR USE:

Clean pacifier before initial use. Single-patient use only.

**Premie:** For 30-40 week gestational age infants and/or for small gestational age babies.

**Newborn:** For babies without teeth (not suggested for the stage of teething).

#### To Clean the CuddleStar:

Hand wash. Air dry. Repeated washing may weaken the fabric and seams. Check frequently and discard if you find any signs of wear. Replace after 3 months of use.

#### To Clean the JollyPop Pacifier :

Clean pacifier before initial use. Open Velcro® closure to remove pacifier. Place pacifier in warm water using a mild liquid soap. To sterilize, place pacifier into boiled water for 2 - 4 minutes. Do not place the CuddleStar in the dishwasher. Replace after 3 months of use.



**Do not tie pacifier around baby's neck as it presents a danger of strangulation.**

3-1502-004C

SandBox Medical LLC

750 Corporate Park, Pembroke, MA 02359 USA  
Tel 877.396.1551 | 781.8266905 Fax 781.826.8210  
orders@sandboxmedical.com sandboxmedical.com

## 7 x 7 Inch Star-Shaped Pacifier Holder with Removable Premie or Newborn Size JollyPop® Pacifier

### PRODUCT SPECIFICATIONS:

Phthalate-free, latex-free, BPA-free and DEHP-free. Manufactured in the USA.

**CuddleStar:** Polyester non-woven outer and inner layer. Velcro® closure for pacifier. New, polypropylene beads contained in inner layer. Meets ASTM F963-11 safety code requirements and is baby safe.

**JollyPop Pacifier:** 100% medical grade silicone. Complies with USA Standard 16CFR 1511.

### INSTRUCTIONS FOR USE:

Clean pacifier before initial use. Single-patient use only.

**Premie:** For 30-40 week gestational age infants and/or for small gestational age babies.

**Newborn:** For babies without teeth (not suggested for the stage of teething).

#### To Clean the CuddleStar:

Hand wash. Air dry. Repeated washing may weaken the fabric and seams. Check frequently and discard if you find any signs of wear. Replace after 3 months of use.

#### To Clean the JollyPop Pacifier :

Clean pacifier before initial use. Open Velcro® closure to remove pacifier. Place pacifier in warm water using a mild liquid soap. To sterilize, place pacifier into boiled water for 2 - 4 minutes. Do not place the CuddleStar in the dishwasher. Replace after 3 months of use.



**Do not tie pacifier around baby's neck as it presents a danger of strangulation.**

3-1502-004C

SandBox Medical LLC

750 Corporate Park, Pembroke, MA 02359 USA  
Tel 877.396.1551 | 781.8266905 Fax 781.826.8210  
orders@sandboxmedical.com sandboxmedical.com