

SandBox Medical

Jellies™

Gel Pillow & Positioning Aid



Meeting All Developmental Needs

6" x 6" Pillows

- Thin: for micropreemies & preemies
- Medium: for fullterm babies
- Thick: for children & adults

6" x 9" Cushions

- Where a larger support is required
- Ideal mattress for micropreemies

New UltraGel!

- Higher viscosity 6" x 9" cushions

Single patient use

Disposable, baby-soft covers available for each size

Made in the USA.



Safe, Affordable, Long-Lasting Gel Support

Give every baby the benefit of gel pillows.

Soft, pliable support for baby's head and body.

- Helps to minimize head molding and skin breakdown.
- Helps to relieve pressure and encourage proper positioning.
- Provides uniform support for the soft cranium of the pre-term infant.
- Nontoxic gel is guaranteed mold and bacteria free. Preservatives in the gel eliminate risk of mold and fungus.
- Lasts 4 to 6 weeks under normal conditions.
- Also makes a great hot or cold pack for nursing moms!

Soft, pliable Jellies help reduce the need for frequent repositioning.



Bright blue color won't get lost in the bed.

For more support and for bigger babies, push the gel to one side, fold under and tape.

Research has shown that gel mattresses are especially effective in reducing vibration during transport.*

* Journal of Perinatology (2015) 35, 110–114; doi:10.1038/jp.2014.172; published online 27 November 2014



Contact us at: orders@sandboxmedical.com

With Jellies™, Every Baby Can Benefit from Gel Support That's Affordable & Long-Lasting.

Gel pillows are an effective way to minimize pressure and provide uniform support for the soft cranium of the pre-term infant. The soft surface and supple support of Jellies Gel Pillows keeps babies more comfortable and minimizes the re-positioning needed to alleviate pressure points.

Easy To Use, Find and Store.

With its bright blue color, Jellies won't get lost in the bed, so you'll save on replacement pillows. You can fold Jellies to adjust the size and thickness. So one size can be used as baby grows. That helps to reduce inventory requirements, ordering and stocking.

Product Specifications

Latex free.

Outer Material: soft polyolefin film.

Contents: Water based non-toxic food grade gel. Methylparaben anti-bacterial agent. See SDS for more information.

Instructions For Use

Single Patient Use Only.

Prewarm gel in incubator or by kneading under warm water (in plastic bag). DO NOT microwave. Cover with a soft cloth or disposable cover. See Jellies complete *Instructions For Use* packaged with product.

Warnings

Jellies are indicated as single patient use to prevent possible cross contamination between patients. Any reuse of Jellies is not warranted by the manufacturer.

Hospital use only. Jellies should only be used in a monitored environment.

DO NOT allow infant to lay face down on pillow. Danger of suffocation. DO NOT X-Ray through pillow. DO NOT microwave.

Jellies comes in thicknesses and sizes to meet all developmental needs.

| Item # | Description | Qty. Per Box |
|------------|--|--------------|
| JE1504212 | 6 x 6 x 3/8 inch Thin Gel Pillow | 12 |
| JE1504312 | 6 x 6 x 5/8 inch Medium Gel Pillow | 12 |
| JE1504512 | 6 x 6 x 3/4 inch Thick Gel Pillow | 12 |
| JE1504412 | 6 x 9 x 1/2 inch Gel Cushion | 12 |
| JE1504812 | 6 x 9 x 3/4 inch Gel Cushion | 12 |
| JE1504835 | 6 x 9 x 3/4 inch Gel Cushion w/ printed cover | 35 |
| JE15041130 | ULTRAGel 6 x 9 x 1 inch Cushion, w/ printed cover | 30 |
| JE15047100 | Disposable printed covers (for 6 x 6 inch) | 100 |
| JE15046100 | Disposable printed covers (for 6 x 9 inch) | 100 |

| Philips Gel-E Donut * | Compares to Jellies: |
|----------------------------|----------------------|
| 92025-A Gel-E Donut Premie | Thin |
| 92025-B Gel-E Donut Small | Medium |
| 92025-C Gel-E Donut | Thick |
| 91033-2 Squishion 2 | 6 x 9 inch |
| 91033-3 Squishion 3 | 6 x 9 inch |

*Gel-E Donut is a registered trademark of Philips Medical.

6 x 6 inch square comes in three thicknesses – to meet the needs of micro-preemies, infants, children and adults.



6 x 9 inch rectangular Jellies – perfect mattress for micro-preemies or hip-roll for bigger babies.



5-1504-006J