




Choose Premie or Newborn Pacifier Size

- 7 x 7 inch star-shaped pacifier holder with removable Premie or Newborn size JollyPop Pacifier
- For hospital and home use
- Single patient use



Made in the USA.



Affordable & Convenient

CuddleStar safely keeps the pacifier close to baby

Helps even tiny babies to keep a pacifier in place.

- Velcro® closure makes it easy to remove and replace JollyPop Pacifier.
- Flexible so you can easily and gently stabilize the pacifier.
- Multiple arms for baby to grasp or cuddle.
- Soft non-woven fabric – appropriate for hospital/home use.
- Partially filled with polypropylene beans – safely contained in a separate inner sack.
- Carefully calculated weight and size keeps the pacifier near baby’s mouth.
- Comes with removable Premie or Newborn size JollyPop® Pacifier (medical grade silicone, latex free).



Soft and pliable CuddleStar is easy for baby to grasp and cuddle.

Tuck the star around baby for side-lying comfort.

Item #	Description	Quantity
VPCS1502112	CuddleStar with JollyPop Premie Pacifier	Box of 12
VPCS1502212	CuddleStar with JollyPop Newborn Pacifier	Box of 12

Product Specifications

Phthalate-free, latex-free, BPA-free and DEHP-free. Manufactured in the USA.

CuddleStar: Polyester non-woven outer and inner layer. Velcro® closure for pacifier.

New, polypropylene beads contained in inner layer. Meets ASTM F963-11 safety code requirements and is baby safe.

JollyPop Pacifier: 100% medical grade silicone. Complies with USA Standard 16CFR 1511.

Instructions For Use

Clean pacifier before initial use. Single-patient use only. See complete Instructions For Use provided with product.

Premie: For 30-40 week gestational age infants and/or for small gestational age babies.

Newborn: For babies without teeth (not suggested for the stage of teething).

 **Do not tie pacifier around baby’s neck as it presents a danger of strangulation.**

NP5-1502-004H

Contact Us For Samples

vizientTM

Awarded Supplier

2014 Pediatric Award Winner

novaplusTM

Novaplus is a registered trademark of Vizient, Inc.

Contract number MS3961